


RUTH LEE™



# FULL BODY CPR TRAINING MANIKIN

## FOR PROFESSIONALS WHO SAVE LIVES

Take CPR training out of the classroom and into the environments where emergencies actually happen. The Ruth Lee Full Body CPR Manikin is a full-size (1.8 m) adult manikin built for realistic, scenario-based training in confined spaces, outdoor settings, and complex rescue scenes.

Supplied with or without a Simulaids Brad torso for accurate CPR training, it allows teams to practise life-saving skills under real-world conditions. Tough, hygienic, and safe to handle, with no pinch points and easy-clean materials, this manikin bridges the gap between textbook training and true operational readiness.

**PROTECT  
MEDICAL**

PART OF PROTECT MEDICAL.  
GET IN TOUCH TO FIND OUT  
MORE ABOUT THE RANGE  
OF MEDICAL AND RESCUE  
PRODUCTS AVAILABLE FROM  
OUR SISTER BRANDS.

SPEAK TO THE EXPERTS:

☎ +44 (0) 1490 413 282 ✉ MAIL@RUTHLEE.CO.UK 💻 WWW.RUTHLEE.COM



# FULL BODY CPR TRAINING MANIKIN

## IDEAL FOR THE FOLLOWING TRAINING SCENARIOS:

First responder Cardiopulmonary Resuscitation (CPR), Confined Space, Road Traffic Collision (RTC), entrapment under a heavy object, or indeed any scenario where CPR could be required under challenging circumstances.

## SPECIAL MATERIALS:

Constructed from the same flame retardant material used in Police ballistic/stab vests, offering 4x the strength and abrasion resistance of 16oz canvas or P.V.C. (conforms to OEKO-TEX® Standard 100, BS EN ISO 12947-2 & 14116: 2008). High density webbing (25mm rated to 200 Kg, 50mm rated to 435Kg with a +/-5% range).

## FEATURES & BENEFITS:

- Full bodied manikin promotes more realistic and effective training.
- Our Full Body CPR manikin comes complete with a Simulaids Brad torso. However, some other torso may fit. Please contact us on 01490 413 282 to check before purchasing if you wish to use a CPR torso you already have .
- Soft body joints eliminate the risk of pinch welts and bruising; common injuries when using plastic bodied manikins.

- Supplied in a tough carrying holdall for easy handling and storage.
- Comes complete with boots and overalls, which are easily replaced if damaged.

## NOT SUITED FOR:

- Use in hot fire training - See our Fire House manikin.
- Use in water – See our Water Rescue range.

## OPTIONAL EXTRAS

- Tough drag protectors and overalls to greatly extend the life of your manikin.

## EXTENSIVELY USED BY:

- UK Emergency Services.
- A wide spectrum of companies across industries.
- Health & Safety and First Aid Training organisations.
- Equivalent organisations across North & South America, Europe, Australia and Asia.

## AVAILABLE WEIGHTS & SIZES:

DESCRIPTION	MODEL CODE	HEIGHT (M)	WEIGHT (KG)	NATO STOCK #
20 Kg Adult with Torso	RLNCPR/TOR	1.8	20	N/A
20 Kg Adult without Torso	RLNCPR	1.8	13	N/A
50 Kg Adult with Torso	RLNCPR/TOR50	1.8	50	N/A
50 Kg Adult without Torso	RLNCPR50	1.8	43	N/A



Please be aware that the simulated "skin" on Simulaids torsos is porous. As a result, contact with dyed or heavily pigmented materials may cause permanent staining. Synthetic blood can also stain the surface. While staining does not affect the functionality of the manikin, it may be considered unsightly over time. Replacement skins can be purchased directly from Simulaids. Ink, dye, and blood stains are not covered under warranty by Simulaids or Ruth Lee Ltd.

Information correct at time of publication 19.12.25 Please be aware material specification is subject to change at any time, sometimes without notice. If specification is vital please contact our Head Office 01490 413 282